

BUB3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BUB3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BUB3L; hBUB3;

Catalog #: A6536

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9184

Isotype: IgG

Swiss Prot: O43684

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein involved in spindle checkpoint function. The encoded protein contains four WD repeat domains and has sequence similarity with the yeast BUB3 protein. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

To place an order, please [Click HERE](#).