

CALU

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CALU

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6538

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 813

Isotype: IgG

Swiss Prot: O43852

Purity: Affinity purification

For research use only.

Background:

The product of this gene is a calcium-binding protein localized in the endoplasmic reticulum (ER) and it is involved in such ER functions as protein folding and sorting. This protein belongs to a family of multiple EF-hand proteins (CERC) that include reticulocalbin, ERC-55, and Cab45 and the product of this gene. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).