

CENPC

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:106kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CENPC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MIF2; hcp-4; CENP-C; CENPC1;

Background:

Centromere protein C 1 is a centromere autoantigen and a component of the inner kinetochore plate. The protein is required for maintaining proper kinetochore size and a timely transition to anaphase. A putative pseudogene exists on chromosome 12.

To place an order, please [Click HERE](#).

Catalog #:A6557

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1060

Isotype:IgG

Swiss Prot:Q03188

Purity:Affinity purification

For research use only.