

CYP4V2

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:61kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CYP4V2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BCD; CYP4AH1;

Background:

This gene encodes a member of the cytochrome P450 hemethiolate protein superfamily which are involved in oxidizing various substrates in the metabolic pathway. It is implicated in the metabolism of fatty acid precursors into n-3 polyunsaturated fatty acids. Mutations in this gene result in Bietti crystalline corneoretinal dystrophy.

To place an order, please [Click HERE](#).

Catalog #:A6573

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:285440

Isotype:IgG

Swiss Prot:Q6ZWL3

Purity:Affinity purification

For research use only.