

DHDH

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DHDH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

2DD; HUM2DD;

Catalog #:A6577

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27294

Isotype:IgG

Swiss Prot:Q9UQ10

Purity:Affinity purification

For research use only.

Background:

This gene encodes an enzyme that belongs to the family of dihydrodiol dehydrogenases, which exist in multiple forms in mammalian tissues and are involved in the metabolism of xenobiotics and sugars. These enzymes catalyze the NADP1-linked oxidation of transdihydrodiols of aromatic hydrocarbons to corresponding catechols. This enzyme is a dimeric dihydrodiol dehydrogenase, and it differs from monomeric dihydrodiol dehydrogenases in its high substrate specificity for trans-dihydrodiols of aromatic hydrocarbons in the oxidative direction.

To place an order, please [Click HERE](#).