

EEF1B2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EEF1B2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EF1B; EEF1B; EEF1B1;

Catalog #: A6580

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1933

Isotype: IgG

Swiss Prot: P24534

Purity: Affinity purification

For research use only.

Background:

This gene encodes a translation elongation factor. The protein is a guanine nucleotide exchange factor involved in the transfer of aminoacylated tRNAs to the ribosome. Alternative splicing results in three transcript variants which differ only in the 5' UTR.

To place an order, please [Click HERE](#).