

EIF2S3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:51kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EIF2S3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EIF2; EIF2G; eIF-2gA; EIF2gamma;

Background:

The protein encoded by this gene is the largest subunit of a heterotrimeric GTP-binding protein involved in the recruitment of methionyl-tRNA(i) to the 40 S ribosomal subunit.

To place an order, please [Click HERE](#).

Catalog #:A6581

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1968

Isotype:IgG

Swiss Prot:P41091

Purity:Affinity purification

For research use only.