

## EIF3I

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**37kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human EIF3I

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

TRIP1; EIF3S2; TRIP-1; PRO2242; eIF3-p36; eIF3-beta;

*To place an order, please [Click HERE](#).*

**Catalog #:**A6582

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8668

**Isotype:**IgG

**Swiss Prot:**Q13347

**Purity:**Affinity purification

For research use only.