

## FCN1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**35kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human FCN1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FCNM;

**Catalog #:**A6587

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2219

**Isotype:**IgG

**Swiss Prot:**O00602

**Purity:**Affinity purification

For research use only.

**Background:**

The ficolin family of proteins are characterized by the presence of a leader peptide, a short N-terminal segment, followed by a collagen-like region, and a C-terminal fibrinogen-like domain. The collagen-like and the fibrinogen-like domains are also found separately in other proteins such as complement protein C1q, C-type lectins known as collectins, and tenascins. However, all these proteins recognize different targets, and are functionally distinct. Ficolin 1 encoded by FCN1 is predominantly expressed in the peripheral blood leukocytes, and has been postulated to function as a plasma protein with elastin-binding activity.

*To place an order, please [Click HERE](#).*