

FGF14

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FGF14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FHF4; FHF-4; SCA27; FGF-14;

Background:

The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. A mutation in this gene is associated with autosomal dominant cerebral ataxia. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A6588

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2259

Isotype:IgG

Swiss Prot:Q92915

Purity:Affinity purification

For research use only.