

FLOT2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:47kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FLOT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ESA; ECS1; ESA1; ECS-1; M17S1;

Background:

Caveolae are small domains on the inner cell membrane involved in vesicular trafficking and signal transduction. This gene encodes a caveolae-associated, integral membrane protein, which is thought to function in neuronal signaling.

To place an order, please [Click HERE](#).

Catalog #:A6590

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2319

Isotype:IgG

Swiss Prot:Q14254

Purity:Affinity purification

For research use only.