

FRZB

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FRZB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FRE; OS1; FZRB; hFIZ; FRITZ; FRP-3; FRZB1; SFRP3; SRFP3; FRZB-1; FRZB-PEN;

Background:

The protein encoded by this gene is a secreted protein that is involved in the regulation of bone development. Defects in this gene are a cause of female-specific osteoarthritis (OA) susceptibility.

To place an order, please [Click HERE](#).

Catalog #:A6592

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2487

Isotype:IgG

Swiss Prot:Q92765

Purity:Affinity purification

For research use only.