

GORASP1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 46kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GORASP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P65; GOLPH5; GRASP65;

Catalog #: A6609

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64689

Isotype: IgG

Swiss Prot: Q9BQQ3

Purity: Affinity purification

For research use only.

Background:

The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is a membrane protein involved in establishing the stacked structure of the Golgi apparatus. It is a caspase-3 substrate, and cleavage of this encoded protein contributes to Golgi fragmentation in apoptosis. This encoded protein can form a complex with the Golgi matrix protein GOLGA2, and this complex binds to the vesicle docking protein p115. Alternative splicing results in multiple transcript variants of this gene.

To place an order, please [Click HERE](#).