

HAS3

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:62kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HAS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6617

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3038

Isotype:IgG

Swiss Prot:O00219

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).