

HYAL2

Reactivity:Human Mouse Rat

Tested applications:WB IF FC

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:200 FC 1:20 - 1:50

Calculated MW:54kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HYAL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LUCA2;

Catalog #:A6624

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8692

Isotype:IgG

Swiss Prot:Q12891

Purity:Affinity purification

For research use only.

Background:

This gene encodes a weak acid-active hyaluronidase. The encoded protein is similar in structure to other more active hyaluronidases. Hyaluronidases degrade hyaluronan, one of the major glycosaminoglycans of the extracellular matrix. Hyaluronan and fragments of hyaluronan are thought to be involved in cell proliferation, migration and differentiation. Although it was previously thought to be a lysosomal hyaluronidase that is active at a pH below 4, the encoded protein is likely a GPI-anchored cell surface protein. This hyaluronidase serves as a receptor for the oncogenic virus Jaagsiekte sheep retrovirus. The gene is one of several related genes in a region of chromosome 3p21.3 associated with tumor suppression. This gene encodes two alternatively spliced transcript variants which differ only in the 5' UTR.

To place an order, please [Click HERE](#).