

IL20

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IL20

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IL-20; IL10D; ZCYTO10;

Catalog #: A6629

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 50604

Isotype: IgG

Swiss Prot: Q9NYY1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis.

To place an order, please [Click HERE](#).