

INPP5J

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:107kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human INPP5J

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PIPP; INPP5; PIB5PA;

To place an order, please [Click HERE](#).

Catalog #:A6631

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27124

Isotype:IgG

Swiss Prot:Q15735

Purity:Affinity purification

For research use only.