

INTS5

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:107kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human INTS5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

INT5; KIAA1698;

Background:

INTS5 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005)

To place an order, please [Click HERE](#).

Catalog #:A6633

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80789

Isotype:IgG

Swiss Prot:Q6P9B9

Purity:Affinity purification

For research use only.