

## IRAK2

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:200

**Calculated MW:**69kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human IRAK2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

IRAK-2;

**Catalog #:**A6635

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3656

**Isotype:**IgG

**Swiss Prot:**O43187

**Purity:**Affinity purification

For research use only.

**Background:**

IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.

*To place an order, please [Click HERE](#).*