

KHDRBS3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KHDRBS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Etle; SALP; SLM2; SLM-2; TSTAR; T-STAR; etoile;

To place an order, please [Click HERE](#).

Catalog #:A6637

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10656

Isotype:IgG

Swiss Prot:O75525

Purity:Affinity purification

For research use only.