

KIF1B

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:204kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KIF1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KLP; CMT2; CMT2A; CMT2A1; HMSNII; NBLST1;

Background:

This gene encodes a motor protein that transports mitochondria and synaptic vesicle precursors. Mutations in this gene cause Charcot-Marie-Tooth disease, type 2A1.

To place an order, please [Click HERE](#).

Catalog #:A6638

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23095

Isotype:IgG

Swiss Prot:O60333

Purity:Affinity purification

For research use only.