

MFGE8

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MFGE8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BA46; HMFG; MFGM; SED1; hP47; EDIL1; MFG-E8; SPAG10; OAcGD3S; HsT19888;

To place an order, please [Click HERE](#).

Catalog #:A6652

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4240

Isotype:IgG

Swiss Prot:Q08431

Purity:Affinity purification

For research use only.