

MYOG

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:25kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human MYOG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MYF4; myf-4; bHLHc3;

Background:

Myogenin is a muscle-specific transcription factor that can induce myogenesis in a variety of cell types in tissue culture. It is a member of a large family of proteins related by sequence homology, the helix-loop-helix (HLH) proteins. It is essential for the development of functional skeletal muscle.

To place an order, please [Click HERE](#).

Catalog #:A6664

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4656

Isotype:IgG

Swiss Prot:P15173

Purity:Affinity purification

For research use only.