

## NHLRC1

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**42kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NHLRC1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

EPM2A; EPM2B; MALIN; bA204B7.2;

**Background:**

The protein encoded by this gene is a single subunit E3 ubiquitin ligase. Laforin is polyubiquitinated by the encoded protein. Defects in this intronless gene lead to an accumulation of laforin and onset of Lafora disease, also known as progressive myoclonic epilepsy type 2 (EPM2).

*To place an order, please [Click HERE](#).*

**Catalog #:**A6669

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**378884

**Isotype:**IgG

**Swiss Prot:**Q6VVB1

**Purity:**Affinity purification

For research use only.