

## NNMT

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**29kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NNMT

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A6673

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4837

**Isotype:**IgG

**Swiss Prot:**P40261

**Purity:**Affinity purification

For research use only.

**Background:**

N-methylation is one method by which drug and other xenobiotic compounds are metabolized by the liver. This gene encodes the protein responsible for this enzymatic activity which uses S-adenosyl methionine as the methyl donor.

**To place an order, please [Click HERE](#).**