

PNPLA3

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PNPLA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADPN; C22orf20; iPLA(2)epsilon;

Catalog #:A6694

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80339

Isotype:IgG

Swiss Prot:Q9NST1

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a triacylglycerol lipase that mediates triacylglycerol hydrolysis in adipocytes. The encoded protein, which appears to be membrane bound, may be involved in the balance of energy usage/storage in adipocytes.

To place an order, please [Click HERE](#).