

POU2AF1

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human POU2AF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BOB1; OBF1; OCAB; OBF-1;

To place an order, please [Click HERE](#).

Catalog #:A6696

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5450

Isotype:IgG

Swiss Prot:Q16633

Purity:Affinity purification

For research use only.