

PRPF19

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:55kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PRPF19

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

el

Synonym:

PSO4; SNEV; PRP19; UBOX4; hPSO4; NMP200;

Background:

PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008)

To place an order, please [Click HERE](#).

Catalog #:A6706

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27339

Isotype:IgG

Swiss Prot:Q9UMS4

Purity:Affinity purification

For research use only.