

QPCT

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:40kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human QPCT

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

QC; GCT; sQC;

Catalog #:A6711

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:25797

Isotype:IgG

Swiss Prot:Q16769

Purity:Affinity purification

For research use only.

Background:

This gene encodes human pituitary glutaminy cyclase, which is responsible for the presence of pyroglutamyl residues in many neuroendocrine peptides. The amino acid sequence of this enzyme is 86% identical to that of bovine glutaminy cyclase.

To place an order, please [Click HERE](#).