

RALB

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RALB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6714

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5899

Isotype: IgG

Swiss Prot: P11234

Purity: Affinity purification

For research use only.

Background:

This gene encodes a GTP-binding protein that belongs to the small GTPase superfamily and Ras family of proteins. GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors.

To place an order, please [Click HERE](#).