

RBP2

Reactivity:Human Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:15kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human RBP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRBP2; RBPC2; CRBP2; CRABP-II;

Catalog #:A6717

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5948

Isotype:IgG

Swiss Prot:P50120

Purity:Affinity purification

For research use only.

Background:

RBP2 is an abundant protein present in the small intestinal epithelium. It is thought to participate in the uptake and/or intracellular metabolism of vitamin A. Vitamin A is a fat-soluble vitamin necessary for growth, reproduction, differentiation of epithelial tissues, and vision. RBP2 may also modulate the supply of retinoic acid to the nuclei of endometrial cells during the menstrual cycle.

To place an order, please [Click HERE](#).