

SARS

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:58kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SARS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SERS; SERRS;

Catalog #:A6733

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6301

Isotype:IgG

Swiss Prot:P49591

Purity:Affinity purification

For research use only.

Background:

This gene belongs to the class II amino-acyl tRNA family. The encoded enzyme catalyzes the transfer of L-serine to tRNA (Ser) and is related to bacterial and yeast counterparts. Multiple alternatively spliced transcript variants have been described but the biological validity of all variants is unknown.

To place an order, please [Click HERE](#).