

STRN3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200 IF 1:50 - 1:100

Calculated MW:87kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human STRN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SG2NA;

To place an order, please [Click HERE](#).

Catalog #:A6756

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:29966

Isotype:IgG

Swiss Prot:Q13033

Purity:Affinity purification

For research use only.