

## TAGLN2

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:** 22kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human TAGLN2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HA1756;

**Catalog #:** A6761

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 8407

**Isotype:** IgG

**Swiss Prot:** P37802

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is similar to the protein transgelin, which is one of the earliest markers of differentiated smooth muscle. The specific function of this protein has not yet been determined, although it is thought to be a tumor suppressor. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*