

TAGLN2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:22kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TAGLN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HA1756;

Catalog #:A6761

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8407

Isotype:IgG

Swiss Prot:P37802

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to the protein transgelin, which is one of the earliest markers of differentiated smooth muscle. The specific function of this protein has not yet been determined, although it is thought to be a tumor suppressor. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).