

TREX1

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:38kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TREX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRV; AGS1; DRN3; HERN3;

Catalog #:A6778

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11277

Isotype:IgG

Swiss Prot:Q9NSU2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a nuclear protein with 3' exonuclease activity. The encoded protein may play a role in DNA repair and serve as a proofreading function for DNA polymerase. Mutations in this gene result in Aicardi-Goutieres syndrome, chilblain lupus, Cree encephalitis, and other diseases of the immune system. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).