

## TXN2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:**18kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human TXN2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MTRX; TRX2; MT-TRX;

**Catalog #:**A6782

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**25828

**Isotype:**IgG

**Swiss Prot:**Q99757

**Purity:**Affinity purification

For research use only.

**Background:**

This nuclear gene encodes a mitochondrial member of the thioredoxin family, a group of small multifunctional redox-active proteins. The encoded protein may play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis.

*To place an order, please [Click HERE](#).*