

WWP2

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 99kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WWP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1

Synonym:

AIP2; WWp2-like;

Catalog #: A6789

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11060

Isotype: IgG

Swiss Prot: O00308

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Nedd4 family of E3 ligases, which play an important role in protein ubiquitination. The encoded protein contains four WW domains and may play a role in multiple processes including chondrogenesis and the regulation of oncogenic signaling pathways via interactions with Smad proteins and the tumor suppressor PTEN. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 10.

To place an order, please [Click HERE](#).