

ZBED1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:78kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ZBED1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALTE; DREF; TRAMP; hDREF;

Background:

This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. It was earlier identified as a gene with similarity to Ac transposable elements, however, was found not to have transposase activity. Later studies show that this gene product is localized in the nucleus and functions as a transcription factor. It binds to DNA elements found in the promoter regions of several genes related to cell proliferation, such as histone H1, hence may have a role in regulating genes related to cell proliferation. Alternatively spliced transcript variants with different 5' untranslated region have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A6792

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9189

Isotype:IgG

Swiss Prot:O96006

Purity:Affinity purification

For research use only.