

ZDHHC17

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:73kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ZDHHC17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HIP3; HYPH; HIP14; HSPC294;

To place an order, please [Click HERE](#).

Catalog #:A6793

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23390

Isotype:IgG

Swiss Prot:Q8IUH5

Purity:Affinity purification

For research use only.