

Bcl9

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:149kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human Bcl9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Gm130; 2610202E01Rik; 8030475K17Rik; A330041G23Rik

To place an order, please [Click HERE](#).

Catalog #:A6795

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:607

Isotype:IgG

Swiss Prot:Q9D219

Purity:Affinity purification

For research use only.