

YBX1

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:1000 IF 1:20 - 1:50

Calculated MW: 36kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human YBX1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

jk

Synonym:

YB1; BP-8; CSDB; DBPB; YB-1; CSDA2; NSEP1; NSEP-1; MDR-NF1

To place an order, please [Click HERE](#).

Catalog #: A6799

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4904

Isotype: IgG

Swiss Prot: P67809

Purity: Affinity purification

For research use only.