

CCS

Reactivity:Human

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:50 - 1:500

Calculated MW:29kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CCS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6803

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9973

Isotype:IgG

Swiss Prot:O14618

Purity:Affinity purification

For research use only.

Background:

Copper chaperone for superoxide dismutase specifically delivers Cu to copper/zinc superoxide dismutase and may activate copper/zinc superoxide dismutase through direct insertion of the Cu cofactor.

To place an order, please [Click HERE](#).