

## BRIP1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**141kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human BRIP1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

OF; BACH1; FANCI;

**Background:**

The protein encoded by this gene is a member of the RecQ DEAH helicase family and interacts with the BRCT repeats of breast cancer, type 1 (BRCA1). The bound complex is important in the normal double-strand break repair function of breast cancer, type 1 (BRCA1). This gene may be a target of germline cancer-inducing mutations.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6804

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**83990

**Isotype:**IgG

**Swiss Prot:**Q9BX63

**Purity:**Affinity purification

For research use only.