

GNL2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:84kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human GNL2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NGP1; Ngp-1; HUMAUANTIG

To place an order, please [Click HERE](#).

Catalog #:A6805

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:29889

Isotype:IgG

Swiss Prot:Q13823

Purity:Affinity purification

For research use only.