

## NTHL1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**34kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NTHL1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FAP3; NTH1; OCTS3; hNTH1;

**Catalog #:**A6820

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4913

**Isotype:**IgG

**Swiss Prot:**P78549

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a DNA N-glycosylase of the endonuclease III family. Like a similar protein in E. coli, the encoded protein has DNA glycosylase activity on DNA substrates containing oxidized pyrimidine residues and has apurinic/aprimidinic lyase activity.

*To place an order, please [Click HERE](#).*