

## PIWIL4

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**97kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PIWIL4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HIWI2; MIWI2;

**Background:**

PIWIL4 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells (Sasaki et al., 2003)

*To place an order, please [Click HERE](#).*

**Catalog #:**A6824

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**143689

**Isotype:**IgG

**Swiss Prot:**Q7Z3Z4

**Purity:**Affinity purification

For research use only.