

RBX1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW:12kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human RBX1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

el

Synonym:

ROC1; RNF75; BA554C12.1;

Catalog #:A6844

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9978

Isotype:IgG

Swiss Prot:P62877

Purity:Affinity purification

For research use only.

Background:

This locus encodes a RING finger-like domain-containing protein. The encoded protein interacts with cullin proteins and likely plays a role in ubiquitination processes necessary for cell cycle progression. This protein may also affect protein turnover. Related pseudogenes exist on chromosomes 2 and 5.

To place an order, please [Click HERE](#).