

## RBX1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:**12kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human RBX1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

el

**Synonym:**

ROC1; RNF75; BA554C12.1;

**Catalog #:**A6844

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9978

**Isotype:**IgG

**Swiss Prot:**P62877

**Purity:**Affinity purification

For research use only.

**Background:**

This locus encodes a RING finger-like domain-containing protein. The encoded protein interacts with cullin proteins and likely plays a role in ubiquitination processes necessary for cell cycle progression. This protein may also affect protein turnover. Related pseudogenes exist on chromosomes 2 and 5.

*To place an order, please [Click HERE](#).*