

WDR48

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WDR48

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P80; UAF1; SPG60;

To place an order, please [Click HERE](#).

Catalog #: A6854

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57599

Isotype: IgG

Swiss Prot: Q8TAF3

Purity: Affinity purification

For research use only.